

A large pile of small, silver-colored metal gears and components, likely made of aluminum or steel. The gears are of various sizes and are scattered across the frame, creating a dense, textured background. The lighting is bright, highlighting the metallic surfaces and the intricate details of the gear teeth.


ACCESSORIES, PARTS AND COMPONENTS

ART 113

ORDER NR

THUMB SCREW ASSEMBLY

for screw type mounting plate or AXIOSPLIT® in the articulator upper or lower member thumb screw assembly 25 mm, head diameter with screw thread, hex 2,0 mm, polyamid and metal washers

**ART 135**


ORDER NR

THUMB SCREW ASSEMBLY FOR FACEBOW**ART 136**

ORDER NR

THUMB SCREW ASSEMBLY FOR BENNETT GUIDANCE


(2 pieces) inside thread 4 mm and outside diameter 9 mm for condylar housings SAM® 2 and SAM® 3

**ART 150**

ORDER NR

BENNETT GUIDANCE, WHITE

for SAM® 2C and SAM® 2P, white (1 Pair), straight guidance, for protrusion and for articulator adjustment

**ART 155**

ORDER NR

REFERENCE POINTER FOR SAM1/SMG**ART 181**

ORDER NR

PROTRUSION INSERT WHITE

1 mm protrusion (2 pieces) for SAM® 2C and SAM® 2P, plastic

**ART 182**

ORDER NR

PROTRUSION INSERT, YELLOW

2 mm protrusion (2 pieces) for SAM® 2C and SAM® 2P, plastic

**ART 183**

ORDER NR

PROTRUSION INSERT, GREEN


3 mm protrusion (2 pieces) for SAM® 2C and SAM® 2P, plastic


**ART 184**

ORDER NR

PROTRUSION INSERT, RED

4 mm protrusion (2 pieces) for SAM® 2C and SAM® 2P, plastic




**PROTRUSION INSERTS, BLUE**

5 mm protrusion (2 pieces) for SAM® 2C and SAM® 2P, plastic

ART 185


ORDER NR

**PROTRUSION INSERTS, BLACK**

6 mm protrusion (2 pieces) for SAM® 2C and SAM® 2P, plastic

ART 186


ORDER NR

**GREASE LUBRICANT**

for metal (50 g) Molykote DX white

ART 192


ORDER NR

**THUMB SCREW FOR CONDYLAR HOUSING**

SAM® 2C and SAM® 2P

ART 220


ORDER NR

**CENTRIC LOCK THUMB SCREW**

for SAM® 2C and SAM® 2P thumb screw with spring loaded contact assembly for condylar housings

ART 245


ORDER NR

**INCISAL PIN IC**

diameter 6 mm x 77 mm long

ART 250


ORDER NR

**THUMB SCREW FOR INCISAL PIN**

for tightening incisal pin, thread 4 mm x 10 mm long, head diameter 12 mm

ART 255


ORDER NR

**THUMB NUT FOR INCISAL TABLE**

with rubber insert in head, thread 4 mm, head diameter 20 mm

ART 265

ORDER NR

**THUMB SCREW FOR INCISAL PIN SLEEVE**

with rubber insert in head, thread 4 mm x 8,5 mm long, head diameter 20 mm

ART 270


ORDER NR

ART 277

ORDER NR

RUBBER INSERT


grey (8 pieces)

**ART 280**

ORDER NR

OCCLUSAL PLANE INDICATOR I C


for SAM® 2C with flat acrylic template insert and angulation scale 0° - 45°, vertical and sagittal adjustment with millimeter scale

**ART 290**

ORDER NR

HEXAGON ANGLE ALLEN WRENCHES


1,5 mm and 2,5 mm

**ART 294**

ORDER NR

ANTERIOR TEETH EVALUATOR I C


for SAM® 2C for upper anterior teeth impression related to axis orbital plane

**ART 305**

ORDER NR

INCISAL PIN IP


for SAM® 2P, diameter 6 mm x 92 mm long

**ART 370**

ORDER NR

INCISAL PIN, SQUARE


for ART 364 / ART 380 / ART 394, non-adjustable in the vertical direction

**ART 435**

ORDER NR

CONDYLAR HOUSING COVERS FOR SAM® 2PX


(1 pair) left and right with thumb screws ART 436, used to hold upper and lower member of the articulator together during eccentric movements.

**ART 436**

ORDER NR

THUMB SCREW FOR CONDYLAR HOUSING COVERS


SAM® 2PX (2 pieces) thread 4 x 7 mm, diameter 9 mm


**ART 521**

ORDER NR

INCISAL PIN SLEEVE III

for SAM® 2P, SAM® 2PX and SAM® 3 includes thumb screws for incisal pin and sleeve



**INCISAL PIN III P**

for SAM® 2P, SAM® 2PX and SAM® 3, with numeric scale, 6 mm diameter x 99 mm long

ART 523

ORDER NR

**END CAPS**

grey, (6 pieces) plastic, for 6 mm round rod

ART 527


ORDER NR

**THUMB SCREW FOR CONDYLAR HOUSING**

SAM® 3 thread 4 x 16 mm, head diameter 16 mm

ART 528


ORDER NR

**TRIPOD ARM**

with end cap grey for articulator upper member

ART 529

ORDER NR

**CONDYLAR HOUSING COVERS FOR SAM® 3**

(1 pair) left and right with thumb screws ART 537, used to hold upper and lower member of the articulator together during eccentric movements.

ART 535


ORDER NR

**THUMB SCREW FOR CONDYLAR HOUSING COVERS**

for SAM® 3 (2 pieces), thread 4 x 7 mm, head diameter 9 mm

ART 537


ORDER NR

**BENNETT GUIDANCE INSERT STRAIGHT**

for SAM® 3, metal (1 pair)

ART 539

ORDER NR

**OCCLUSAL PLANE INSERT, CURVED**

replacement acrylic template, engraved, also used for complete denture setup.

ART 580

ORDER NR

**OPEN-END WRENCH**

size 6 mm, drop-forged, wrench angle 15°

ART 583


ORDER NR

ART 584

ORDER NR

OPEN END WRENCH


size 10, drop-forged, angle of the wrench 15°

**ASP 452**

ORDER NR

AXIOSPLIT® RUBBER RING


for ASP 450 / ASP 451 (2 pieces), white, inside diameter 58 mm and height 21 mm

**ASP 475**

ORDER NR

AXIOSPLIT® NORM INCISAL PIN P


for SAM® 2P and SAM® 3, to set factory incisal pin vertical dimension value in an AXIOSPLIT® equipped articulator, for use with ASP 470

**ASP 480**

ORDER NR

AXIOSPLIT® MAGNET

magnets (2 pieces) with protecting ring for bayonet remover

**ASP 530120**

ORDER NR

AXIOSPLIT® - MAGNET S WITH SOCKET


for AXIOSPLIT®-S plate

**ATB110**

ORDER NR

HEX WRENCH FOR ATB AND AXO


for 7 mm diameter hex screw

**ATB 120**

ORDER NR

THUMB NUT

anodized

**ATB 121**


ORDER NR

ENCLOSURE PLATE**ATB 125**

ORDER NR

SCREW

slot headed, for horizontal rod of the facebow transfer fork assembly, 4 x 6 mm






EARPIECES,

(1 pair) 15 mm, black, version for the Transfer Bow AXIOQUICK® I and II

ATB 201

ORDER NR




WEDGE SHAPED SLIDING BLOCK

black

ATB 230002

ORDER NR




VERTICAL ROD

replacement for ATB 205 and ATB 305 without horizontal rod, clamp or transfer fork, vertical rod length is 68 mm

ATB 306

ORDER NR



HORIZONTAL ROD

replacement for ATB 100, ATB 205 or ATB 305

ATB 307

ORDER NR




EARPIECES

for ATB (2 pieces) black, 25 mm, version for Transfer Bow AXIOQUICK® III

ATB 311


ORDER NR



REPLACEMENT PARTS TRANSFER STAND AX

ATB 398010

ORDER NR



RULER FOR AXIOGRAPH®

steel, 200 mm used to measure AXIOGRAPH® interflag distance

AXO 105

ORDER NR




AXIOGRAPH® FRENCH CURVE

plastic, for manual drawing of Bennett curves from AXIOGRAPH® data output

AXO 106

ORDER NR



AXIOGRAPH® CONVERSION TABLE

used to calculate the Bennett angle, analog curvature, and excursion width for selection of Bennett guidance analog inserts

AXO109


ORDER NR

AXO 130

ORDER NR

CROSS CLAMP


adjustable

**AXO 159**

ORDER NR

RUBBER BAND


grey for AXIOGRAPH® (2 pieces), for the fixation on the neck

**AXO 180**

ORDER NR

MEASUREMENT LOUPE


for AXIOGRAPH®, used to determine condylar analog inclination and pathway shape

**AXO 201**

ORDER NR

CROSSROD


for AXIOGRAPH® recording bow, steel, 6.0 mm diameter, hollow 275 mm long

**AXO 265**

ORDER NR

SIDEARM SLEEVE, BLACK, WITH AXIS PIN


axis pin AXO 280 has a pointed tip and a groove in the shaft

**AXO 270**

ORDER NR

SIDEARM SLEEVE, SILVER, WITH AXIS PIN


axis pin AXO 280 has a pointed tip and a groove in the shaft

**AXO 275**

ORDER NR

SIDEARM SLEEVE, RED, WITH AXIS PIN


axis pin AXO 281 has a flat tip and no groove in the shaft

**AXO 276**

ORDER NR

THUMB SCREW FOR SLEEVE, BLACK


for AXO 265, thread 3 x 6.5 mm, head diameter 7 mm

**AXO 277**

ORDER NR

THUMB SCREW FOR SLEEVE, SILVER


for AXO 270 thread 3 x 6.5 mm, head diameter 7 mm






THUMB SCREW FOR SLEEVE, RED
for AXO 275 thread 3 x 6,5 mm, head diameter 7 mm

AXO 278
ORDER NR




AXIS PIN, POINTED TIP, WITH GROOVE
for AXO 265 / AXO 270 for black or silver AXIOGRAPH® sidearm sleeve

AXO 280
ORDER NR




AXIS PIN, FLAT TIP
for AXO 275 for red AXIOGRAPH® sidearm sleeve

AXO 281
ORDER NR




MEASURING LOUPE FOR AXIOGRAPH®
with writing device with lead 0,5 mm, for measuring Bennett sideshift and recording mandibular movement pathways

AXO 290
ORDER NR




SIDEARM SLEEVE, YELLOW, WITH WRITING DEVICE
with lead holder (0,5 mm), for AXIOGRAPH® lower sidearms

AXO 310
ORDER NR




BENNETT ANALOG INSERT SELECTOR
with angle selection, orange, overlay, with angle indication, for selecting Bennett guidance analog insert from AXIOGRAM graphic output

AXO 315
ORDER NR




NT-CLAMP UPGRADE KIT
for AXIOGRAPH® conversion from double clam to NT Clamp

AXO 330
ORDER NR



NT-CLAMP AND CROSSROD
for AXIOGRAPH®

AXO 331
ORDER NR



PARAOCCLUSAL CLUTCH
AXIODAPT® with hex screw driver


AXO 333
ORDER NR

AXO 341

ORDER NR

HEX SCREW DRIVER


for AXIODAPT®, AXO 333 / 334, 2,0 mm diameter

**AXO 410001**

ORDER NR

SIDEARM LEFT


lower for AXIOGRAPH® III, adjustable horizontally and vertically

**AXO 410002**


ORDER NR

SIDEARM RIGHT

lower for AXIOGRAPH® III, adjustable horizontally and vertically

**AXO 410003**


ORDER NR

NT- CLAMP**AXO 410005**

ORDER NR

NT TOGGLE NUT


complete

**AXO 410006**

ORDER NR

WING SCREW

with moveable toggle

**AXO 410008**

ORDER NR

O-RING

for toggle

**AXO 410009**

ORDER NR

NEW SILVER AXIS


for NT-clamp


**AXO 420**

ORDER NR

SIDEARM WITH FRONT ROD

for AXIOGRAPH® III upper flag bow, left and right sidearms, front crossrod, and recording flag screws






NASION RELATOR
for AXIOGRAPH® III adjustable

AXO 425

ORDER NR


PARTS



VERTICAL HEAD SUPPORT
for AXIOGRAPH® III

AXO 430

ORDER NR



EARPIECE ALIGNMENT FLAG, LEFT
for AXIOGRAPH® III with ear piece and axis pin for axis alignment tube

AXO 440


ORDER NR



EARPIECE ALIGNMENT FLAG, RIGHT
for AXIOGRAPH® III with ear piece and axis pin for axis alignment tube

AXO 441


ORDER NR



RECORDING FLAG, LEFT
for AXIOGRAPH® III

AXO 442


ORDER NR



RECORDING FLAG, RIGHT
for AXIOGRAPH® III

AXO 443


ORDER NR



THUMB SCREW FOR UPPER SIDEARMS
AXIOGRAPH® III (2 pieces) thread 4 x 8 mm,
head diameter 9 mm

AXO 444


ORDER NR



HINGE AXIS TRANSFER CLAMP
for use with upper recording flags to make a hinge axis face-bow transfer to the hinge axis mounting device using the inside surface of the recording flags

AXO 450

ORDER NR



TRANSFER FORK ASSEMBLY II - NT
for AXIOGRAPH® III, with NT-Clamp and transfer fork, vertical rod length is 95 mm

AXO 460


ORDER NR

AXO 470

ORDER NR

ORBITAL REFERENCE PLANE INDICATOR


for AXIOGRAPH® III

**AXO 475**

ORDER NR

ORBITAL REFERENCE POINTER


for AXIOGRAPH® III

**AXO 480**

ORDER NR

VERTICAL ROD WITH MASTOID SUPPORT


for AXIOGRAPH® III

**AXO 485**

ORDER NR

AXIS ALIGNMENT TUBE


for AXIOGRAPH® III (2 pieces)

**AXO 487**

ORDER NR

HORIZONTAL EXTENSION RODS


AXIOGRAPH® III (2 pieces) used to extend upper flagbow side-arms for alignment and stabilization of upper flag bow

**AXO 490**

ORDER NR

ALIGNMENT GAUGE


for AXIOGRAPH® III, used to align the upper flagbow parallel to midsagittal plane and collinear with the hinge axis

**AXO 512**


ORDER NR

THUMB SCREW FOR LOWER SIDEARMS

for thumb screw and sidearm sleeve and dial gauge, thread screw 3 mm x 10 mm, thumb nut 7 mm diameter

**AXO 550**


ORDER NR

RR-CPS KIT FOR AXIOGRAPH III**MOH 100**

ORDER NR

MOUNTING SUPPORT INSTRUMENT

with transfer fork support acrylic (ATB 335)



MOH 109

ORDER NR

THUMB SCREW ASSEMBLY FOR MOH

thumb nut and plastic screw assembly for mounting support instrument




MOH 120

ORDER NR

MOUNTING CAST HOLDER

for centric mounting with 30 flexible plastic rods




MOH 250

ORDER NR

LAB MOUNTING DEVICE FOR MOUNTING SUPPORT INSTRUMENT

for upper cast mounting in a mounting support instrument without using SAM® 2 Articulator upper member, for use with screw type mounting plate or AXIOSPLIT® System.




MOH 265M

ORDER NR

LAB MOUNTING DEVICE FOR MOUNTING SUPPORT INSTRUMENT, MPS

for upper cast mounting in a mounting support instrument without using SAM® 2 Articulator upper member, for use with magnetic mounting plates




MOH 300K

ORDER NR

MOUNTING SUPPORT INSTRUMENT LAB KIT

includes a mounting support instrument, telescopic transfer fork support, centric mounting cast holder with 30 flexible plastic rods




MOH 400K

ORDER NR

MOUNTING SUPPORT INSTRUMENT LAB KIT, PROFESSIONAL

includes a mounting support instrument, telescopic transfer fork support, a centric mounting cast holder with 30 flexible plastic rods and a lab mounting device MOH 250 for use with screw type mounting plate or AXIOSPLIT® System




MOH 415K

ORDER NR

MOUNTING SUPPORT INSTRUMENT LAB KIT, PROFESSIONAL, MPS

includes a mounting support instrument, telescopic transfer fork support, a centric mounting cast holder with 30 flexible plastic rods and a lab mounting device MOH 265, for use with screw type mounting plate or AXIOSPLIT® System




MOH 562

ORDER NR

MOUNTING CAST HOLDER

for mounting stand MOH 562, 30 flexible plastic rods



MPI 205

ORDER NR

MANDIBULAR POSITION INDICATOR CALIBRATION TOOL FOR DIAL GAUGE

used to check and set dial gauge to ZERO, set distance between measuring cube and medial wall




MPI 350

ORDER NR

INCISAL PIN
universal, for MPI**MPS 111**


ORDER NR

MPS ALIGNMENT TOOL
used to centralize the MPS Magnet Assembly on the upper or lower member of the articulator**MRI 311**


ORDER NR

DRILL FOR MRI II
diameter 3 mm, 2,33 mm shaft for use in straight dental hand-piece**MSF 201**


ORDER NR

MAGNETIC SPLITCAST FORMER PLATE, STANDARD, RED
used for making gypsum split casts**MSF 203**


ORDER NR

RUBBER RING FOR MSF STANDARD / REVERSE, BLACK
(2 pieces) inside diameter 79 mm and height 21 mm**MSF 204**


ORDER NR

RUBBER RING STANDARD, RED
(2 pieces) extra height with inside diameter 79 mm and height 24 mm**MSF 251**


ORDER NR


MAGNETIC SPLITCAST FORMER PLATE, REVERSE BLACK
used for making the secondary gypsum split casts prior to pouring the impression. When used in a similar fashion as the regular splitcast, the primary splitcast is automatically made at the time the cast is poured, saving one step.**MSF 301**

ORDER NR

MAGNETIC SPLITCAST FORMER PLATE, SMALL, YELLOW
used for making gypsum split casts**MSF 303**

ORDER NR


RUBBER RING FOR MSF SMALL, WHITE
(2 pieces) inside diameter 68 mm and height 21 mm

**BOXING AID DISTANCE MEASURING DEVICE**

metal gauge for determining minimum gypsum cast thickness in a boxing aid device

MSF 410


ORDER NR

**MAGNETIC SPLITCAST FORMER PLATE, XL**

green, extra large, for making gypsum split casts for full dentures procedures

MSF 451


ORDER NR

**BOXING AID RING, XL**

white plastic, includes green clamp

MSF 452


ORDER NR

**RUBBER RING FOR MSF XL, BLACK**

(2 pieces) inside diameter 102 mm and height 21 mm

MSF 453


ORDER NR

**DOWEL PIN DRILL BIT**

1,6 mm (3 pieces) for drilling dowel pin holes in gypsum casts to receive conical or AXIOPIN parallel dowel pins

MSF 503

ORDER NR

**DOWEL PIN DRILL DEPTH GAUGE FOR AXIODRILL®**

used to establish correct pin drilling depth in gypsum casts

MSF 606

ORDER NR