


## Stops & Hooks



# Surgical and Rectangular Arch Hooks

- Prevent failure of bracket attachments to teeth
- 2 mm cast with hook
- Split hooks can be placed without removing arch wire
- Provide an efficient means of ligation
- Use sliding hook crimping plier for strongest crimp
- Surgical hooks may be used for intermaxillary fixation following surgery
- 10 pieces per package



Surgical hook



Rectangular arch hook

## Surgical Hooks

## Rectangular Arch Hooks

		Surgical Hooks		Rectangular Arch Hooks		
		Slide on	Split	Slide on	Split Right	Split Left
Inches	mm					
.018 x .025	0.46 x 0.64	851-620C	851-622C	851-632C	851-630CR	851-630CL
.021 x .025	0.53 x 0.64	851-621C	851-623C	851-633C	851-631CR	851-631CL

## Sliding Hook Crimping Plier

- Two grooved inserts maintain a firm grip on the tube
- Used with American surgical hooks and rectangular arch hooks

For surgical hooks (shown) 100-107  
 For arch hooks 100-107A



Sliding hook crimping plier

# Stops and Hooks for Round Wire

## Crimpable


- 1 mm and 2 mm lengths
- Wall thickness .006 inches (.15 mm)
- Ideal for light wires
- Use dull hard wire cutter for placement
- Split crimpable can be placed without removing arch wire

## Weldable

- 2 mm lengths
- Wall thickness .015 inches (.38 mm)
- Welded to arch wire prior to ligation



Arch hook







Round arch stop

		Crimpable			Split Crimpable	Weldable	
		Arch hook 10/pkg	1 mm Stop 100/pkg	2 mm Stop 100/pkg	2 mm Stop 50/pkg	Arch hook 10/pkg	2 mm Stop 100/pkg
Inches	mm						
.020	0.51		851-740X	851-740			
.022	0.56	851-701	851-741X	851-741		851-720	
.025	0.64	851-702	851-742X	851-742			
.030	0.76	851-703	851-743X	851-743	851-634		
.036	0.91	851-704	851-744X	851-744		851-723	851-763
.040	1.02	851-705	851-745X	851-745		851-724	851-764
.045	1.14	851-706	851-746X	851-746		851-725	851-765
.051	1.3	851-707	851-747X	851-747		851-726	851-766

## Stops for Rectangular Wire

- All rectangular stops have a wall thickness of .015 inches (.38 mm)
- Split crimpable can be placed without removing arch wire

		Cast Split Crimpable	Split Crimpable	Weldable	
		2 mm Stop	2 mm Stop	1.25 mm Stop	2 mm Stop
		10/pkg	100/pkg	100/pkg	100/pkg
					
Inches	mm				
.018 x .025	0.46 x 0.64	851-635C	851-612	851-614	851-610
.021 x .028	0.53 x 0.71	851-636C	851-613	851-615	851-611


Rectangular arch stop



## Monkey Hook

- .018 linkable S-shaped hook
- Directs impacted teeth into position
- Assists in the correction of severely rotated teeth
- May be hung from a bondable coil eyelet #292-186
- Can be used in multiples to create a chain
- 10 pieces per package


851-640




Bondable coil eyelet

## Arch Lock

- Removable lock secures to arch wires from .012 (0.30 mm) up to .022 x .028 (0.56 x 0.71 mm) and retainer wire up to .030 (0.76 mm)
- Pre-assembled unit allows placement in one piece



Arch lock without hook



Arch lock with hook



Correct placement of arch lock on arch wire

**KIT** without hooks 851-650

- 10 locks without hooks
- 1 hex key handle
- 2 hex keys

**KIT** with hooks 851-651

- 10 locks with hooks
- 1 hex key handle
- 2 hex keys

## Components

- |                        |         |
|------------------------|---------|
| 10 locks without hooks | 851-652 |
| 10 locks with hooks    | 851-653 |
| 1 hex key handle       | 851-654 |
| 5 hex keys             | 851-655 |